


# Achieve Low Flow Sensitivities with Micro Flow Chromatography on the QTRAP<sup>®</sup> 6500 System for Targeted Quantitative Proteomics

## Using the Eksigent ekspert<sup>™</sup> nanoLC 425 System

Christie Hunter, Kelli Jonakin  
AB SCIEX, USA

There has been a significant amount of research focused on discovering proteins/peptides that are differentially expressed in specific cell and disease conditions. To confirm or refute their ultimate utility, many more samples must be analyzed with increased throughput and robustness, which means faster chromatography and/or higher flow rates. At higher flow rates sensitivity may decrease due to reduced ionization efficiency compared to nanospray ionization. Thus the best sensitivities are typically achieved using nL/min flow rates. There has been increased interest in working in the microflow regime (3-50  $\mu$ L/min) to obtain a good balance between throughput, robustness and sensitivity. Coupled with recent advances in IonDrive<sup>™</sup> technology that provides higher sensitivity detection, microflow chromatography should provide a step forward in productivity and ease of use at a good sensitivity level.



**Figure 1. Comparing Signal Intensities of Standard Protein Digest.** A beta-galactosidase digest was analyzed by nanoflow LC on a QTRAP<sup>®</sup> 5500 System and compared to microflow LC on the QTRAP<sup>®</sup> 6500 System. Similar signal intensities were observed with similar separation quality but with >2x faster total run times with microflow LC.



The purpose of this study is to explore the level of sensitivity and speed of analyses achievable using two hybrid triple quadrupole linear ion trap systems operating at different flow regimes. The goal is to better understand the most efficient strategy for larger scale protein quantification studies.

## Key Features of the Microflow QTRAP<sup>®</sup> 6500 System for Peptide Quantitation

- IonDrive<sup>™</sup> Technology on the QTRAP<sup>®</sup> 6500 system provides:
  - Increased ionization efficiency and heat transfer with the new IonDrive<sup>™</sup> Turbo V source (Figure 2)
  - Increased ion sampling efficiency and ruggedness with the new IonDrive<sup>™</sup> QJet ion guide
  - Increased dynamic range with new IonDrive<sup>™</sup> high energy detector technology
- Mass range of  $m/z$  5 – 2,000 provides versatility for peptide quantitation
- The ekspert<sup>™</sup> nanoLC 425 system has the flexibility to support a broad range of flow rates (from nano to microflow) in a single system
  - Plug and Play flow modules
  - High speed autosampler enables excellent injection reproducibility with little or no sample waste


## Experimental

**Sample Preparation:** Beta-Galactosidase digest and 6 Peptide mixture was obtained from AB SCIEX. The Six Protein Digest was obtained from Michrom BioResources. Crashed plasma matrix was prepared by mixing equal volumes of plasma and acetonitrile, followed by centrifugation.

**Chromatography:** All chromatography was performed on the Eksigent ekspert™ nanoLC 425 System, using the 0.1–1  $\mu\text{L}/\text{min}$  or 1–10  $\mu\text{L}/\text{min}$  flow modules. Nanoflow chromatography was performed using the Eksigent cHiPLC® System in trap elute mode with 75  $\mu\text{m}$  x 15 cm, C18 Eksigent ChromXP™ cHiPLC columns at 300 nL/min flow rate. Microflow chromatography was performed using 300  $\mu\text{m}$  x 15 cm, C18 ChromXP columns at 4  $\mu\text{L}/\text{min}$  flow rate. Injection volumes used ranged from 2 – 5  $\mu\text{L}$ .

**Mass Spectrometry:** NanoLC experiments were performed on the QTRAP® 5500 System with NanoSpray® Source. MicroLC experiments were performed on the QTRAP® 6500 System with IonDrive™ Turbo V Source with 25 $\mu\text{m}$  ID hybrid electrodes. MRM transitions were optimized for each peptide and used on both instruments. Standard concentration curves were performed to evaluate impact of flow rates and separation times on the two different MS systems. All samples were analyzed in triplicate.

**Data Processing:** Lower limits of quantification (LLOQ) were determined using MultiQuant™ Software.



**Figure 3. Signal Intensity Comparison of Eight Tryptic Peptides from Six Protein Mixture at Concentrations near LLOQ.** Signal at 7.6 amol on column of the protein digest is shown for both the QTRAP® 5500 and the QTRAP 6500 System. In this case, the signal intensities found in the microflow case (bottom) showed a small improvement over the nanoflow data (top).



**Figure 2: IonDrive™ Turbo V Source.** The optimized geometry and larger diameter heaters provide more efficient heat transfer to a larger cross section of the spray cone. Equipped with the lower inner diameter hybrid electrodes from Eksigent, this source provides a robust easy to optimize solution for microflow chromatography.

## MicroFlow LC for Peptide Quantification


When moving to higher flow rates on the QTRAP® 6500 System, the IonDrive™ TurboV Source provides high efficiency ionization and increased ruggedness (Figure 2). For best performance at microflow rates, the sources are plumbed with the hybrid electrodes specifically designed for microflow<sup>1</sup>. These electrodes significantly reduce post column dead-volume for reduced dispersion and sharper peak widths. In this work, the 25  $\mu\text{m}$  ID electrode was used, which is ideal for 300  $\mu\text{m}$  ID columns and a flow rate range from 3 – 25  $\mu\text{L}$ .

Previous work explored the difference in sensitivities between different nano and micro flow rates by performing standard concentration curves on a set of tryptic peptides.<sup>2</sup> It was found that there is roughly a 3x loss in sensitivity when moving from nanoflow to the higher 4  $\mu\text{L}/\text{min}$  microflow rate. Other peptide quantitation studies have measured a roughly 3-5 fold increase in sensitivity when moving from the QTRAP® 5500 System to the QTRAP® 6500. From these two studies it was hypothesized that nanoflow rate peptide quantification assays currently performed on the QTRAP 5500 system could be moved to microflow on the QTRAP 6500 system, resulting in similar sensitivities but with improved robustness and increased throughput.

## Comparing the Sensitivity Differences Across LCMS Platforms

The first comparative experiment analyzed the signal intensities for a standard beta-galactosidase digest at 10 fmol on column. The signal intensities were quite comparable between the nanoflow experiment on the QTRAP 5500 system and the microflow experiment on the QTRAP 6500 System (Figure 1). The total run time was reduced by 2 fold in the microflow experiment but separation resolution was preserved.

The next experiment monitored a set of tryptic peptides in a simple matrix using concentrations curves to assess sensitivity. Again, The LLOQs were determined using both LCMS systems and the results across the peptides were compared (Table 1). While some variation is observed across the peptide group, on average a 2x lower LLOQ was observed when using microflow configuration on the QTRAP 6500 system. Figure 3 shows the MRM signal difference across the 8 peptides observed by the two techniques, at an on-column concentration just above the LLOQs. Representative data is shown for one of the peptides from carbonic anhydrase (Figure 4). The LLOQ was found to be at 1.9 amol for the microflow experiment on the QTRAP 6500 system, ~2x more sensitive than that observed on the QTRAP 5500 system at nanoflow (3.8 amol on column).



**Figure 4. LLOQ Comparisons for Peptide DGPLTGTYR from Carbonic Anhydrase.** The MRM signal at the LLOQ for both the nanoflow and microflow shows a ~2x improvement in sensitivity when using the microflow chromatography on the QTRAP® 6500 System.

The final comparison experiment was performed using the Bradykinin peptide in protein precipitated plasma. Standard concentration curves were run on both LCMS platforms and LLOQs determined. The results for the microflow experiment on the QTRAP 6500 system provided an LLOQ of about 6 amol on column (Figure 5). On the QTRAP 5500 system at nanoflow, the LLOQ was found to be 12.5 amol. Again roughly equivalent or slightly better sensitivity was achieved when using microflow on the QTRAP 6500 System.

**Table 1: Lower Limits of Quantification (LLOQ) Obtained for Eight Tryptic Peptides on the two LCMS Systems.** Standard concentration curves in simple matrix were generated and the LLOQs were determined using both nanoflow on QTRAP 5500 system and microflow on QTRAP 6500 system. The results for the peptides show some variation observed across peptides but on average a 2x lower LLOQ was seen on the microflow QTRAP 6500 system.

Peptide	Fragment Ions Summed for Quant	QTRAP® 6500 System – 4 μL/min LLOQ (amol on column)	QTRAP® 5500 System – 300 nL/min LLOQ (amol on column)	Sensitivity Improvement on Microflow QTRAP 6500 System
IDALNENK	2y4, 2y6, 2y7	.48	1.9	4.0
TPEVDDEALEK	2y102+, 2y7, 2y8	3.8	3.8	1.0
VLVLDTDYK	2y5, 2y6, 2y7	1.9	3.8	2.0
AEFVEVTK	2y5, 2y6	3.8	3.8	1.0
ATEEQLK	2y5	7.6	3.8	0.5
DGPLTGTYR	2y5	1.9	3.8	2.0
VGDANPALQK	2y5, 2y6, 2y7	.95	3.8	4.0
VLDALDSIK	2y6, 2y7, 2y8	1.9	3.8	2.0
<b>Average Difference:</b>				<b>2.1</b>

## Conclusions


The QTRAP® 6500 System provides high sensitivity, excellent robustness and broad dynamic range for MRM quantification of peptides. The flexibility and reproducibility of the Eksigent ekspert™ nanoLC 425 makes it an ideal LC system for labs performing a broad range of proteomics workflows, which may include both nano and micro flow rate applications.

Microflow chromatography on the QTRAP® 6500 system provides an easy to use, high-throughput workflow for performing targeted peptide quantitation. Assays currently performed by nanoflow LC on the QTRAP 5500 system could be easily translated to the microflow QTRAP 6500 system for accelerated sample analysis with similar sensitivities.

- Robust LCMS with microflow rates, using both direct injection and trap-elute configurations
- Achieve nanoflow sensitivities with microflow on a QTRAP® 6500 system
- Increase throughput 2X or more
- Fast, high precision autosampler for microflow chromatography

## References

1. Higher Sensitivity and Improved Resolution Microflow UHPLC with Small Diameter Turbo V™ Source Electrodes. AB SCIEX technical note 4590211-01
2. Exploring the Sensitivity Differences for Peptide Quantification in the Low Flow Rate Regime - Eksigent ekspert™ nanoLC 400 System. AB SCIEX technical note 6560212\_02.



Row	Component Name	Actual Conc.	Num. V	Mean	Standard Devi	Percent CV	Accuracy
1	Bradykinin.2y7	6.30	3 of 3	5.411e0	8.429e-1	15.58	85.88
2	Bradykinin.2y7	12.50	3 of 3	1.263e1	1.120e0	8.87	101.01
3	Bradykinin.2y7	25.00	3 of 3	2.555e1	2.747e0	10.75	102.18
4	Bradykinin.2y7	50.00	3 of 3	4.626e1	5.584e0	12.07	92.52
5	Bradykinin.2y7	100.00	3 of 3	1.069e2	6.212e0	5.81	106.83
6	Bradykinin.2y7	200.00	3 of 3	1.942e2	4.001e0	2.06	97.09
7	Bradykinin.2y7	400.00	3 of 3	4.279e2	8.055e0	1.88	106.98
8	Bradykinin.2y7	800.00	3 of 3	8.413e2	9.923e0	1.18	105.16
9	Bradykinin.2y7	4000.00	3 of 3	4.007e3	8.415e1	2.10	100.19
10	Bradykinin.2y7	20000.00	3 of 3	2.056e4	2.363e2	1.15	102.78
11	Bradykinin.2y7	100000.00	3 of 3	9.937e4	1.641e3	1.65	99.37

**Figure 5. Quantification of Bradykinin using Microflow LC on QTRAP® 6500 System.** Standard concentration curve of Bradykinin peptide in protein precipitated plasma provided an LLOQ of 6.3 amol on column. Linearity was very good across the limited dynamic range interrogated (>4 orders in this example). The equivalent experiment using nanoflow chromatography on the QTRAP 5500 system resulted in an LLOQ of 12.5 amol on column for the same y7 fragment ion.

For Research Use Only. Not for use in diagnostic procedures.

© 2013 AB SCIEX. The trademarks mentioned herein are the property of AB Sciex Pte. Ltd. or their respective owners. AB SCIEX™ is being used under license.

Publication number: 7480213-01



**Headquarters**  
500 Old Connecticut Path | Framingham, MA 01701 USA  
Phone 508-383-7700  
[www.absciex.com](http://www.absciex.com)

**International Sales**  
For our office locations please call the division headquarters or refer to our website at [www.absciex.com/offices](http://www.absciex.com/offices)