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SCIEX Triple Quad™ 4500 LC-MS/MS System on Simultaneously Detection of Trimethylamine and Trimethylamine Oxide in human serum

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SCIEX China Clinical BU

Keywords: SCIEX Triple Quad™ 4500; Trimethylamine; Trimethylamine Oxide, Human Specimen

Abstract: Objective: The aim of this work is to develop and validate a rapid, sensitive and accurate method to determine trimethylamine N-oxide (TMAO) and trimethylamine (TMA) in human serum that can reflect the risk of coronary atherosclerotic heart diseases (CAD) using stable isotope dilution liquid chromatography tandem mass spectrometry (LC-MS/MS). Method: The analytes were separated by an ACQUITY UPLC BEH Amide, 2.1 × 100 mm, 1.7 μm and was analyzed via electrospray ionization mass spectrometry. A binary gradient by mixing water and acetonitrile with formic acid and ammonium formate was generated to elute the analytes. Aliquots of 50 μl standard solutions, serum samples and quality controls were mixed with 20 μl internal standards and were extracted with 130 μl acetonitrile: methanol= 85:15. The mix was vortexed and centrifuged (12000 rpm for 10 min) before the supernatant was transferred and analyzed. Result: The chromatographic separation method was controlled within 5 min. The mean linear correlation coefficients were higher than 0.99. The analytical recoveries ranged from 92.0 to 110.9% with an average of 103.4%. The intra-run and total CVs for measurement of these metabolites were 0.84-6.19% and 1.71-8.26% respectively. Conclusion: An isotope dilution LC-MS/MS method for simultaneous measurement of two CAD related metabolites was developed and validated in this work.

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- [1] Senthong V, Wang Z, Fan Y, Wu Y, Hazen SL, Tang WH. Trimethylamine N-Oxide and Mortality Risk in Patients With Peripheral Artery Disease[J]. Journal of the American Heart Association. 2016,5(10).
- [2] Nie J, Xie L, Zhao BX, Li Y, Qiu B, Zhu F, Li GF, He M, Wang Y, Wang B et al. Serum Trimethylamine N-Oxide Concentration Is Positively Associated With First Stroke in Hypertensive Patients[J]. Stroke. 2018,49(9):2021-8.
- [3] Ge X, Zheng L, Zhuang R, Yu P, Xu Z, Liu G, Xi X, Zhou X, Fan H. The Gut Microbial Metabolite Trimethylamine N-Oxide and Hypertension Risk: A Systematic Review and Dose-Response Meta-analysis[J]. Advances in nutrition (Bethesda, Md). 2020,11(1):66-76.

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