

SCIEX Dry Pump MSR 130

OPERATOR GUIDE

Disclaimer

This document is provided to customers who have purchased SCIEX equipment to use in the operation of such SCIEX equipment. This document is copyright protected and any reproduction of this document or any part of this document is strictly prohibited, except as SCIEX may authorize in writing.

Software that may be described in this document is furnished under a license agreement. It is against the law to copy, modify, or distribute the software on any medium, except as specifically allowed in the license agreement. Furthermore, the license agreement may prohibit the software from being disassembled, reverse engineered, or decompiled for any purpose. Warranties are as stated therein.

Portions of this document may make reference to other manufacturers and/or their products, which may contain parts whose names are registered as trademarks and/or function as trademarks of their respective owners. Any such use is intended only to designate those manufacturers' products as supplied by SCIEX for incorporation into its equipment and does not imply any right and/or license to use or permit others to use such manufacturers' and/or their product names as trademarks.

SCIEX warranties are limited to those express warranties provided at the time of sale or license of its products and are the sole and exclusive representations, warranties, and obligations of SCIEX. SCIEX makes no other warranty of any kind whatsoever, expressed or implied, including without limitation, warranties of merchantability or fitness for a particular purpose, whether arising from a statute or otherwise in law or from a course of dealing or usage of trade, all of which are expressly disclaimed, and assumes no responsibility or contingent liability, including indirect or consequential damages, for any use by the purchaser or for any adverse circumstances arising therefrom.

Trademarks and/or registered trademarks mentioned herein, including associated logos, are the property of AB Sciex Pte. Ltd., or their respective owners, in the United States and/or certain other countries (see sciex.com/trademarks).

For Research Use Only. Not for use in Diagnostic Procedures

AB Sciex[™] is being used under license.



Made in Czechia Manufactured for: AB Sciex LLC 500 Old Connecticut Path Framingham, Massachusetts 01701 USA

Copyright notice

© 2023 DH Tech. Dev. Pte. Ltd. Published: 9/4/2023



Contents

1.	Safety and compliance	6
	1.1. Definition of Warnings and Cautions	6
	1.2. Safety symbols	7
2.	General description.	
	2.1. Overview	8
	2.2. Interfaces and control panel	
	2.2.1. Auxiliary connector socket	10
	2.2.2. Logic interface	
	2.2.3. Auto-run	
	2.2.4. Pump controller	12
3	Technical data	. 13
0.	3.1. Operating and storage conditions.	-
	3.1. Operating and storage conditions. 3.2. Performance.	
	3.3. Materials exposed to gases.	
	3.4. Intake gas specification. 3.5. Electrical data.	
		15
4.	Installation.	. 16
	4.1. Installation safety.	16
	4.2. Unpack and inspect.	16
	4.3. Mechanical installation.	17
	4.4. Leak test the system	20
	4.5. Electrical installation.	
	4.5.1. Connect the electrical supply	20
	4.6. Connect the system.	21
	4.6.1. Connection for remote control and monitoring.	21
Б	Commission the nump	. 22
5.		. 22
6.	Operation.	. 23
	6.1. Operation safety	23
	6.2. Operational modes	24
	6.2.1. Manual control mode	24
	6.3. Logic interface data	25
	6.4. Parallel control and monitoring	27
	6.5. Analogue speed control	28
	6.6. Hardware configuration	30
	6.7. Operational modes	30
	6.8. Start the pump	30
	6.9. Shut down the pump	
	6.10. Restart the pump	31
	6.10. Restart the pump. 6.11. Unplanned shutdown.	
		31
	Commission the pump	. 2
6.	•	
	6.2. Operational modes	24
	6.2.1. Manual control mode	24
	6.3. Logic interface data	25
	6.4. Parallel control and monitoring	27
	6.5. Analogue speed control.	28
	6.6. Hardware configuration.	30
	6.7. Operational modes	30
	•	
	6.10. Restart the pump	31
	6.11. Unplanned shutdown	31
	6.11. Unplanned shutdown	31 31



7. Maintenance and service.	
7.1. General maintenance	
7.2. Maintenance plan	34
7.3. Inspect and clean the inlet strainer	34
7.4. Clean the ventilation slots	
7.5. Electrical safety check	34
7.6. Return the equipment or components for service (SCIEX)	35
8. Fault finding	
8.1. Alarm indicator codes	37
9. Storage	38
10. Disposal	39



List of Figures

Figure 1: Pump features	9
Figure 2: Control panel	10
Figure 3: Valve connector.	11
Figure 4: Dimensions	13
Figure 5: Performance graph	14
Figure 6: Unpacking Instruction.	17
Figure 7: Logic interface connections - parallel control	28
Figure 8: Logic interface connections - analogue speed control	29
Figure 9: Analogue speed control	29
Figure 10: Maximum vapour handling capacity	32
Figure 11: Gas ballast connection	33



1. Safety and compliance

For safe operation from the start, read these instructions carefully before you install or commission the equipment and keep them safe for future use. Read all the safety instructions in this section and the rest of this manual carefully and make sure that you obey these instructions.

The instruction manual is an important safety document that we often deliver digitally. It is your responsibility to keep the instruction manual available and visible while working with the equipment. Please download the digital version of the instruction manual for use on your device or print it if a device will not be available.

1.1. Definition of Warnings and Cautions

For safe operation from the start, read these instructions carefully before you install or commission the equipment and keep them safe for future use.

Read all the safety instructions in this section and the rest of this manual carefully and make sure that you obey these instructions. The equipment must only be operated and maintained by trained personnel in the proper condition and as described in this instruction manual.

Obey local and state requirements and regulations. If you have any questions about safety, operation or maintenance of the device, please contact our nearest subsidiary.

Important safety information is highlighted as warning and caution instructions which are defined as follows. Different symbols are used according to the type of hazard.

WARNING:

If you do not obey a warning, there is a risk of injury or death.

CAUTION:

If you do not obey a caution, there is a risk of damage to equipment, related equipment or process.

NOTICE:

Information about properties or instructions for an action which, if ignored, will cause damage to the equipment.

We reserve the right to change the design and the stated data. The illustrations are not binding.



1.2. Safety symbols

The safety symbols on the products show the areas where care and attention is necessary.

The safety symbols that we use on the product or in the product documentation have the following meanings:

Warning/Caution Risk of injury and/or damage to equipment. An appropriate safety instruction must be followed or a potential hazard exists.
Warning - Heavy object Risk of injury or damage to equipment. Identifies a possible hazard from a heavy object.
Warning - Dangerous voltage Risk of injury. Identifies possible sources of hazardous electrical shock.
Warning - Hot surfaces Risk of injury. Identifies a surface capable of inflicting burns through contact.
Symbol - Protective earth Identifies an electrical equipment earth (ground) terminal.
Warning - Use protective equipment Risk of injury. Use appropriate Personal Protective Equipment (PPE) when performing the task.



2. General description

2.1. Overview



WARNING: INCORRECT USE OF EQUIPMENT

Risk of injury or damage to equipment. Incorrect use of the equipment can cause injury. The user is responsible for the safe operation, installation and monitoring of the system.

CAUTION: GAS BALLAST INSTALLATION

Risk of damage to the equipment. If water vapour is present within the system or is pumped by the product the gas ballast must be used to prevent damage to the product. The gas ballast should continue to be applied with the pump running, for at least 1 hour after the water vapour source has been removed to enable the pump to dry.

The MSR 130 is a dry, multi-stage roots vacuum pump that offers high pumping speed in a compact form. The pump is connected by a power cord and is designed for use on clean duty applications.

The cord rating is:

Supply voltage (V)	Current (A)	Frequency (Hz)
200 - 240	16	50 or 60

The pump is not designed for use with flammable, corrosive, toxic or other hazardous gases. Gas or oxygen can mix in the pump system.

The dry pump operates in steady-state gas load conditions with an inlet pressure from 1 - 15 mbar. If the inlet pressure is more than 15 mbar for an extended period, the pump controller reduces the motor speed. The motor speed is increased when the pressure is reduced, or when the transient overload protection has recovered. Refer to *Table: Performance data* for information on maximum inlet pressure.

Refer to *Figure: Pump features*. The system is supplied with NW40 inlet port (item 3) and NW25 exhaust port (item 2).

The system is air cooled by a fan installed in the pump enclosure. The fan continues to operate for 5 minutes after the pump is set to off. The pump has a thermal protection device that will stop the motor if a thermal overload, for example in high ambient temperature, occurs. You must restart the pump after it has cooled down.

The pump mechanism is operated by an electric motor driven by an internal pump controller. The rotational elements of the pumping mechanism are simply supported at the end by lubricated bearings. PFPE lubricant is used and this is contained in both the gearbox and motor ends of the pump module. These are "sealed for life" and the lubricant does not need to be replaced until the pump service interval is reached.

The system has a gas ballast connection, if it is necessary to purge the system.

The pump is designed to remove the residual gases used in high vacuum systems. The gases are:

- Air
- Oxygen (O₂) < 21% by volume
- Neon
- Nitrogen
- Argon




Helium

You can use the pump to pump water vapour. You must use the gas ballast when water vapour is pumped. The water vapour must not condense in the pump.

If you wish to pump a gas not contained within the list above, or if you are unsure about safe gas mixtures or concentration levels, contact your supplier for advice. If you do not contact your supplier, you may invalidate the warranty on the pump. The pump is not suitable for pumping aggressive, corrosive or mixture of flammable gases.

Figure 1 Pump features




- 3. Inlet port
- 5. Leveling feet x 4
- 7. Lifting eye positions
- 9. Mains connector port
- 11. 15-way D-type connector

- 4. Castor x 4
- 6. Gas ballast connection
- 8. Mains circuit breaker
- 10. Protective earth stud
- 12. Auxiliary connector



2.2. Interfaces and control panel

Figure 2 Control panel



- 1. Start button
- 3. Standby increase button
- 5. Service indicator LED
- 7. Auto-run indicator LED
- 9. Standby increase indicator LED
- 11. Stop button
- 13. Power status LED

- 2. Standby button
- 4. Standby decrease button
- 6. Alarm indicator LED
- 8. Standby decrease indicator LED
- 10. Standby indicator LED
- 12. Run indicator LED

2.2.1. Auxiliary connector socket

An auxiliary control connection on the rear panel controls an optional inlet valve. This inlet valve can be operated in parallel with the normal pump output signal. Refer to *Figure: Pump features* item 6. The valve is usually closed and:

- will open when the normal signal is active (pump at speed)
- will close when you select stop, or if there is a fault condition.


The reaction time will be in line with the valve selection. The output signal is 24 V d.c. Refer to *Figure: Valve connector* for the polarity of the connector pins when the connector is energized.

The auxiliary connector is regulated to 24 V d.c. to control the accessories. Refer to *Table: Auxiliary load currents*, the output will shut down to protect the pump controller.

Refer to *Table: Recommended mating plugs* for the recommended mating plugs.



Figure 3 Valve connector



Pin number	Signal	Polarity
1	Valve 1 – Control output	Positive
2	Valve 2 – Control output	Positive
3	Valve 1 - Return	Negative
4	Valve 2 - Return	Negative

Table 1 Auxiliary load currents

Description	Data
Connector plug	Phoenix part number SACC-DSI-M 8FS-4CON-M12/0.5
Voltage output	24 V d.c25%, +10% (18 V d.c. to 26.4 V d.c.)
Output power	2 Channels with 4 watts per channel

Table 2 Recommended mating plugs

Mating connector plug	Phoenix part number
Screw connection, straight	SACC-M 8MS-4CON-M-SW
Solder connection, straight	SACC-M 8MS-4CON-M
Screw connection, right angle	SACC-M 8MR-4CON-M-SW
Solder connection, right angle	SACC-M 8MR-4CON-M

2.2.2. Logic interface

The logic interface is designed to support the serial control and the parallel control to monitor and control through one connector.

The pump controller can be operated through the 15-way D-type logic interface connector (refer to *Figure: Pump features*, item 8). The signals on the logic interface are:

- Control inputs: switch type and analogue signals that control the pump.
- Status outputs: to identify the status of the system.

For serial control select RS232 or RS485. Refer to Serial Comms Interface manual.

For control modes refer to Table: Manual control mode.

For logic interface data refer to Logic interface data.

2.2.3. Auto-run

The auto-run setting configures the pump to start automatically when the power is switched on, without any customer intervention.



You can configure the auto-run through serial communications, or with the start or stop button. Push and hold the start or stop button for more than eight seconds to enable or disable the auto-run setting. The auto-run LED displays the auto-run setting.

The pump can be stopped by either manual, parallel or serial control modes when in auto-run mode.

Note:

Refer to Operation safety for guidelines on how to configure different operational modes and any associated earnings to consider.

2.2.4. Pump controller



CAUTION: ELECTRICAL SUPPLY

Risk of damage to equipment. Do not disconnect the pump from the electrical supply until it has has fully stopped.

The pump controller contains the drive electronics to control the pump operation.

The pump controller controls the supply of electric current to the motor for the operating conditions. This allows the pump to be connected to a single phase mains supply.

The interface control panel is installed on the pump controller. The pump can be operated:

- Manually, with the buttons on the user interface control panel. Refer to Figure: Control panel.
- Remotely, with the serial communications or the digital and analogue process control (parallel). The control is connected to the pump through the 15-way D-type interface connector. Refer to *Figure: Pump features* and *Connection for remote control and monitoring*.

3. Technical data

3.1. Operating and storage conditions

Table 3 Operating and storage conditions

Range	Data
Ambient operating temperature range	5 °C to 40 °C
Ambient operating humidity range	10% to 90% RH (non-condensing)
Maximum operating altitude	2000 m
Ambient storage temperature range	-30 °C to 70 °C
Pollution degree	2
Equipment type	Indoor use
Enclosure protection	IP20

Note:

The pump-system is designed to meet the electrical safety creepage and clearance requirements for 3000 m altitude. For maximum thermal performance, de-rating is required between 2000 m and 3000 m. Reduce the maximum operating ambient temperature from 40 °C to 35 °C between 2000 m and 3000 m linearly. This de-rating meets the requirements of IEC 282-1.2.

Figure 4 Dimensions




Table 4 Mechanical data

Parameter	Units	Data
Overall dimensions (L x W x H)	mm	696 x 297 x 343
Maximum tilt angle	deg	25°
Nominal rotational speed	rpm	9000
Mass	kg	77
Inlet connection	-	NW40
Outlet connection	-	NW25



3.2. Performance


Table 5 Performance data


Parameter	Units	Data
Maximum pumping speed with 16.6 mm restrictor em- bedded	m ³ /h	126
Maximum pumping speed at atmospheric pressure	m ³ /h	6
Ultimate pressure	mbar	0.02
Maximum continuous inlet pressure	mbar	20
Maximum continuous exhaust pressure	bar	0.4
Suck-back protection	-	By exhaust valve
Leak tightness (static)	mbar l/s ⁻¹	< 1 x 10 ⁻⁵
Maximum water vapor pumping rate	kg/h	0.3

* Values for the maximum pumping speed are taken at a sea level. At levels higher than sea level, the pump performance can be affected and the rotational speed can drop below the nominal value.

Table 6 Sound data

Declared dual-number noise-emission values in accordance with ISO 4871	
	Data
Measured A-weighted emission sound pressure level, L _{pA} at ultimate vacuum 1 m from the pump in free space dB(A)	55.4
Uncertainty, K _{pA} dB(A)	2.5
Values from ISO 3744: 2010	







3.3. Materials exposed to gases

The materials and the components exposed to the gases in the pump are:

- Iron
- Aluminium alloy
- Steel
- Fluoroelastomer (seals)
- PTFE (seals)
- PFPE lubricant

3.4. Intake gas specification

Table 7 Intake gas specification

Specification	Reference data
Intake gas	Dry air, nitrogen, argon or other inert gases (except Xenon and Radon)
Intake gas interface	8 mm tube compression fitting
Maximum allowable intake gas supply pres- sure	7 bar (gauge), 100 psig, 7 x 10 ⁵ Pa
Minimum allowable intake gas supply pres- sure	3 bar (gauge), 45 psig, 3 x 10 ⁵ Pa
Minimum purity	ISO 8573 solids class 1
Typical flow at 3 bar	4 slm

Refer to Gas ballast operation on page 31.

3.5. Electrical data

Table 8 Electrical rating for continuous operation

Supply voltage	Phase	Frequency	Rated current	Max continuous	Max overload
(V a.c. rms)		(Hz)	(A rms)	current (A rms)	current (A rms)
200 - 240 +/-10%	Single	50 - 60	8.5	10	16.5 for maximum 3 minutes

Note:

In a transient overload condition the maximum overload current drawn by the pump is 10 amps. Refer to Transient overload protection for more information.

4. Installation

4.1. Installation safety



CAUTION: SAFETY INSTRUCTIONS

Follow all safety instructions and take note of all appropriate precautions.



CAUTION: BLOCKED EXHAUST PIPELINE

Risk of damage to equipment. Make sure the exhaust pipeline is not blocked. If an exhaust isolation valve is used, make sure that the pump is not operated with the valve closed.

Possible hazards on the dry pumping system include electricity, process chemicals, and Fomblin® (PFPE) oil:

- Contact us or see sciex.com/support/request-support for advice or assistance on installation.
- Do not remove the temporary covers from the system inlet and exhaust until ready to connect.
- Do not operate the system unless the inlet and exhaust are connected to the vacuum and exhaust extraction system.
- Isolate the other components in the process system from the electrical supply to prevent accidental operation.
- Electrical supplies are potentially hazardous energy sources. Lockout and tagout before you do the maintenance.
- Obey all national and local rules and safety regulations when you install the system.
- Tighten the cables, hoses and pipework during installation to prevent a trip hazard.
- Make sure that the installation area is clean and free from debris and contamination before you install the pump.
- Make sure that all facilities given in this manual are available for the system to perform correctly.

4.2. Unpack and inspect



WARNING: DAMAGED PUMP

Risk of injury to people or damage to equipment. Do not use the pump if it is damaged. Failure to do so can result in injury to people and/or damage to equipment.




CAUTION: PUMP LIFE

Risk of damage to equipment. Unpack the pump carefully and avoid excessive shock to the pump. Excessive shock will damage the bearings and reduce the life of the pump.



Figure 6 Unpacking Instruction



Refer to Table: Mechanical data for the mass of the pump.

To unpack the pump follow the instruction as shown in the Figure: Unpacking Instruction.

The pump is attached to a wooden pallet with metal brackets in a cardboard package. The metal brackets must be removed and kept.

- 1. Open the cardboard box from the top.
- 2. Remove the packing material. Keep all the packing materials for use in inspection and if the pump is returned for service.
- 3. Examine the pump for damage. If the pump is damaged, tell your supplier and the carrier immediately. Give the supplier and the carrier the information that follows:
 - part number of the pump
 - serial number of the pump
 - order number
 - supplier's invoice number
- 4. Do not use the pump if the pump is damaged.
- 5. If the pump is not to be used immediately, store the pump in the conditions given in Storage.
- 6. The pump is supplied with sealed inlets and outlet to prevent the entry of dust and vapor. Do not remove the seals until the pump is ready to be installed on the vacuum system.

4.3. Mechanical installation



WARNING: HEAVY OBJECT

Risk of physical injury. Use suitable equipment to lift the pump.



WARNING: TRIPPING HAZARD

Risk of injury or damage to equipment. Make sure that cables or pipe work attached to the pump are installed carefully to prevent a slip or trip hazard and the damage to the cable.



WARNING: LIFTING EYE INTEGRITY

Risk of injury or damage to equipment. Make sure that the maximum angle between the paired slings used to lift the system is 45 degrees. Increased angle could decrease the lifting eye integrity.



WARNING: INSTALLATION SAFETY

Risk of injury or damage to equipment. Install the pump in the vacuum system before you connect the pump to the power supply. This will make sure that the pump is not operated and injure people during installation.



WARNING: TOPPLE HAZARD

Risk of injury or damage to equipment. Do not exceed the topple angle when you move the pump or the pump can fall over and injure people.



WARNING: HEAVY OBJECT

Push-pull hazard with strain injury. Push or pull the pump system only for short distance and over flat surfaces. Lift the system if the floor is uneven or has obstacles.



WARNING: EXHAUST BLOCKAGE

Risk of injury or damage to equipment. A peak pressure of 3 bar(g) can be generated in the pump if the exhaust or pipework attached to the exhaust are blocked.



WARNING: STATIC PRESSURE

Risk of injury or damage to equipment. Do not step or stand on the pump, it is not designed to withstand large static loads.



CAUTION: INSTALLATION SAFETY

Risk of damage to the environment. Obey all local legislation when the pump is installed or removed to reduce the impact of the pump on the environment.



CAUTION: CONDENSATE DRAINAGE

Risk of damage to equipment. Use a catchpot to prevent the drainage of condensate back into the system. Condensate that drains back into the system could damage the pump.

Note: Do not touch power cable with hot surfaces.

Obey the instructions that follow when you install the pump:

- Make sure that there is no blockage to access the pump electrical supply cable or the other controls.
- Make sure that there is a minimum air gap of 40 mm on all sides of the pump to allow effective air circulation.



- Make sure that the system is on a firm and levelled surface that can support the mass of the pump.
- Make sure that the system is installed away from combustible materials.
- You must do a risk assessment of the location and make sure that you can move the pump safely and as per the local and national manual handling guidelines.

To attach the system to the floor:

- 1. Put the M10 (class 12.9) bolts (not supplied) in the slotted holes of the shipping brackets to attach the the system to the floor.
- 2. Tighten the M10 (class 12.9) bolts.
- 3. If the vibration transmission to the floor is a problem, install the applicable vibration isolators (not supplied) between the brackets and the tie bolt or the stud.
- 4. Level the pump to a maximum of 5 degrees (measured at the pump inlet).

Note:

The shipping brackets that are used to attach the pump to the shipping pallet can be used for seismic restraint. The shipping brackets are designed to withstand a level 4 earthquake at a ground floor installation.

To move the pump:

- 1. Attach the correct lifting equipment to all 3 lifting eyebolts. Refer to Figure: Pump features.
- 2. Move the pump carefully to the installation area. Refer to *Table: Mechanical data* for weight of the pump.
- 3. Adjust the levelling feet (Refer to *Figure: Pump features*) to make sure that the dry pumping system is levelled and not supported by the castors. The recommended jacking height is 5 mm.
- 4. Remove the plastic caps from the inlet and exhaust before you connect the pump to the vacuum system. Use appropriate NW vacuum fittings for the connection to the system.

Obey the instructions that follow when you connect the pump to the vacuum system:

- Connect the pump to an exhaust line to minimize the noise and the exhaust emissions.
- Make sure that the pipeline connected to the pump inlet is as short as possible. Make sure that the pipeline has a minimum internal diameter to get maximum pump speeds.
- Put support under the vacuum pipeline to prevent the load on the coupling joints.
- Make sure that the pump exhaust line is not blocked as a pressure of 3 bar(g) can be generated in the exhaust pipework. Connect the pump with appropriate pipework and fittings.
- If necessary, install flexible bellows in the system pipelines to reduce the transmission of the vibration and to prevent the load on the coupling joints. The pressure rating of the bellows must be higher than the highest pressure generated in the system. We recommend that you use the manufacturer's bellows.
- Install an inlet isolation valve in the pipeline between the vacuum system and the pump. The inlet isolation valve isolates the vacuum system from the pump when the system is set to off. This prevent the flow of the process gases and particles back into the vacuum system.
- Make sure that the sealing surfaces are clean and are not scratched.
- We recommend that you use an exhaust extraction system appropriate for use with all pumped process gases. Make sure that the exhaust extraction system is not blocked or obstructed when the pump is in operation.



4.4. Leak test the system



WARNING: SYSTEM LEAK TEST

Risk of injury or damage to equipment. Do the leak test of the system after installation. Seal all the leaks found to prevent leakage of dangerous substances out of the system and leakage of air into the system.

We will accept no liability or warranty claims for damages caused from flammable mixtures because of air leaks.

- 1. Do the leak test of the vacuum system after installation.
- 2. Seal all the leaks found.

4.5. Electrical installation



WARNING: ELECTRICAL CONNECTION

Risk of electric shock. The electrical installation must be done by a qualified person. Always make the electrical connections to the pump after the pump has been installed on the vacuum system.



WARNING: INSTALLATION SAFETY

Risk of electric shock or damage to equipment. The pump must be electrically installed in accordance with regional and local codes, and must obey the local and national safety requirements.



WARNING: HAZARDOUS VOLTAGE

Risk of electric shock or damage to equipment. The logic interface is 30 V maximum rated PELV and must only be connected to PELV interfaces. Failure to use an correctly rated supply could result in electric shock.

The pump is a cord-connected device rated at 8.5 A and must be installed in accordance with local electrical regulations. The pump is to be supplied from a single phase 200 - 240 V branch circuit protected supply rated at 15 A maximum. Refer to *Table: Recommended mating plugs*.

4.5.1. Connect the electrical supply



WARNING: PROTECTIVE EARTH CONNECTION

Risk of electric shock. Make sure that the pump and electrical cables are suitably protected against earth (ground) faults. We recommend you attach a protective earth (ground) conductor (with a cross sectional area of 2.5 mm²/14 AWG) to the protective earth (ground) stud.

Make the electrical connection to the pump mains connector port (*Figure: Pump features*) with an appropriate cord set. Refer to *Table: Recommended mating plugs*.

The pump must be grounded through the earth (ground) conductor of the IEC60320 connector.



4.6. Connect the system

4.6.1. Connection for remote control and monitoring

To operate the pump with parallel or serial control, use the 15-way D-type connector. Refer to *Figure: Pump features*.

Refer to *Figure: Logic interface connections - parallel control, Figure: Logic interface connections - analogue speed control* for details of the logic interface pins.



5. Commission the pump

To commission the pump:

- 1. Make sure that all openings to atmospheric pressure in the foreline vacuum system are closed.
- 2. Connect the power cord.
- Set the mains circuit breaker (*Figure: Pump features*, (8)) to ON. Make sure that the power indicator LED (*Figure: Control panel*, (14)) illuminates. If the LED does not illuminate, contact us.
- 4. Press the start button (*Figure: Control panel*, (1)) until the run indicator LED (*Figure: Control panel*, (13)) begins to flash.
- 5. When the dry pumping system starts and continues to operate, if an alarm condition is indicated:
 - a. Shut down the dry pumping system. Refer to *Shut down the pump* on page 30.
 - b. Refer to *Fault finding* on page 36.
 - c. If the problem has not been rectified, contact us.
- 6. After you commission the dry pumping system:
 - a. to continue to operate the system, refer to Start the pump on page 30
 - b. to shut down the system, refer to *Shut down the pump* on page 30.



6. Operation

6.1. Operation safety



WARNING: DAMAGED PARTS

Risk of electric shock. Do not operate the pump with any parts of the enclosures removed or damaged as there can be a risk of an electric shock.



WARNING: OPERATIONAL SAFETY

Risk of injury or death. Do not expose any part of the human body to vacuum as this can result in injury or death of people.



WARNING: PUMP TEMPERATURE

Risk of injury. Do not touch the pump inlet manifold or exhaust when the pump is in operation as the pump inlet manifold and exhaust can get hot. Allow the pump to cool down after the pump has stopped.



WARNING: AUTOMATIC START

Risk of injury or damage to equipment. The system has an auto-run mode which, if configured, is designed to automatically start the pump system once power is applied.



CAUTION: CONDENSATE DEPOSITION

Risk of damage to equipment. Do not use the pump to pump particulates or condensate. Deposition may occur in the pump. This can degrade the pump performance and reduce the pump life.



CAUTION: MAXIMUM CONTINUOUS OPERATING PRESSURE

Risk of damage to equipment. The pump can be controlled, stopped or put in standby mode only through the interface used to start the pump. In auto-run mode, the stop button on the user interface panel will override the starting interface and stop the pump.



CAUTION: EN55011

This equipment is not intended for use in residential environments and may not provide adequate protection to radio reception in such environments.

Note:

This is an industrial (Class A, Group 1) product as defined by EN55011 and EN61325. To ensure compliance with European Electromagnetic Compatibility (EMC) requirements for EMC emissions, please note that it is not intended for use in domestic buildings or in properties directly connected to an electrical supply network which also supplies domestic buildings. "Group 1" is defined as equipment which does not use RF energy as an intrinsic part of operation or process.



6.2. Operational modes

The pump can be controlled through:

- Manual control mode with the buttons on the user interface panel
- Parallel control mode through the 15-way D-type logic interface connector located at the rear of the pump
- Serial control mode through the 15-way D-type logic interface connector located at the rear of the pump
- Auto-run

6.2.1. Manual control mode

Refer to Figure: Control panel for the pump control functions.

Table 9 Manual control mode

Operation	Button	Status
Apply power	Mains power	The pump will remain off (factory default). The power indicator will illuminate.
Start the pump	Start button	The pump will accelerate up to full run- ning speed. The run indicator will flash while the pump accelerates. The run indicator will remain on when the pump reaches full speed.
Select and deselect the standby speed	Standby mode select button	When engaged, the standby indicator will illuminate and the pump will run at the standby speed setting. The pump is set by default at 70% of full speed.
Increase or decrease the pump speed when in standby mode	Standby speed in- crease button	The pump speed will increase. The increase standby indicator will re- main illuminated when the pump reaches a maximum of 100% of full speed.
	Standby speed decrease button	The pump speed will decrease. The decrease standby indicator will re- main illuminated when the pump reaches a minimum of 50% of full speed.
Select and deselect the Auto-run function	Start or stop button (>8 sec)	When engaged, the auto-run indicator will illuminate. The pump will restart auto- matically after the power has been re- stored.

Start and stop

Use the buttons as shown in *Figure: Control panel* to start and stop the pump.

Note:

The stop command does not isolate the pump from the electrical supply.



Standby

In standby mode the pump operates at a reduced speed to improve the service life of the pump.

- 1. Push the standby button to select standby mode. The pump will run at factory default standby speed (87% of full speed).
- 2. Adjust the speed with the increase and decrease standby speed buttons. The maximum standby speed is 100% of the default run speed and the minimum standby speed is 83% of the default run speed.
- 3. A single short push will change the speed by 1% of the default run speed. Hold the button to change the speed by 1% per second.
- 4. Once adjusted, the pump will return to the new user defined speed each time standby speed is selected.
- 5. Push the standby button to return to normal run speed.

Table 10 LED Indicators

Description	Function
Power indicator	Indicates that electrical mains supply to the pump is ON.
Run indicator	Indicates that the pump is running.
LED continuously ON	Pump runs at full speed
LED flashing	When changing speed
LED OFF	Pump is not running
Standby mode indicator	Indicates that the standby mode has been selected.
Standby speed increase in- dicator	The indicator will blink with every short push of the standby speed increase button. The indicator will remain ON when maximum standby speed has been reached.
Standby speed decrease indi- cator	The indicator will blink with every short push of the standby speed decrease button. The indicator will remain ON when minimum standby speed has been reached.
Auto-run indicator	Indicates that the auto-run mode has been selected.
Standby indicator LED	Once the standby function is selected the run indicator will flash at 50% duty while the rotational speed is changed. Once steady state running speed has been reached both the run speed and standby indicators will be illuminated.
Service indicator	Indicates that a service interval has been reached.
Alarm indicator	Indicates an alarm has been triggered.

Refer to Figure: Control panel for LED locations on the pump dashboard.

6.3. Logic interface data

The pumps have a 15-way D-type logic interface connector located on the user interface panel (*Figure: Pump features*). The logic interface connector can be plugged directly into the 200 W Turbo Instrument Controller (TIC) with a pump TIC cable and in conjunction with controller extension cables.

For Turbo controller, or Turbo and Active Gauge controller (TAG), a suitable connector mating half must be used (not supplied) to connect the pump to the customer control system. Refer to *Table: Logic interface pins* for the electrical connections.



Table 11 Interface technical data

Connector	15-way D-type (male)			
Start, serial enable and remote enable:				
Enable control voltage: low (closed) Disable control voltage: high (open)	0 to 0.8 V d.c. (IOUT = 0.55 mA nominal) 4 to 26.4 V d.c. (internal pull-up to a Thevenin equivalent circuit: 5.3 V and 11 kOhms nomi- nal)			
Standby control input:				
Enable control voltage: low (closed) Disable control voltage: high (open)	0 to 0.8 V d.c. (IOUT = 0.3 mA nominal) 4 to 26.4 V d.c. (internal pull-up to a Thevenin equivalent circuit: 2.5 V and 10.3 kOhms nominal)			
Analogue and RS485 enable control input	S:			
Enable control voltage: low (closed) Disable control voltage: high (open)	0 to 0.8 V d.c. (IOUT = 0.55 mA nominal) 4 to 26.4 V d.c. (internal pull-up to a Thevenin equivalent circuit: 5.5 V and 11.4 kOhms nominal)			
Analogue speed input	8.4 to 10 V d.c. directly proportional to the mo- tor speed e.g. 8.4 V = 125 Hz, 10 V = 150 Hz			
Voltage accuracy	± 5% full scale			
NORMAL status output:				
Type < Normal speed (default 80%) > Normal speed Maximum current rating Maximum voltage rating	Open collector transistor plus pull up resistor OFF (4.7 k pull up + diode to 12 V d.c.) ON (< 0.8 V d.c. sinking 10 mA) 10 mA 28.8 V d.c			
FAIL status output:				
Analogue 10 V reference	+ 10 V d.c. analogue voltage reference uni- polar output with diode protection			
Voltage accuracy output	± 2% full scale			
Current	± 5 mA for specified accuracy			

Table 12 Logic interface pins

Pin num- ber	Signal	Polarity	Use
1	Analogue speed enable- control input	-	Connect to pin 2 (0 V) to enable analogue speed control via pin 9.
2	0 V control reference	-	0 V reference for all control and status sig- nals listed within this table.
3	Start/Stop – control input	-	Connect to pin 2 (0 V) to start the pump system.



Pin num- ber	Signal	Polarity	Use
4	Standby – control input/ serial-RX/RS-485 A-	-	Connect to pin 2 (0 V) to enable standby speed when the serial enable control input is inactive.
5	Serial enable – control input	-	Connect to pin 2 (0 V) to enable serial com- munications.
6	RS-232/RS-485 – control input	-	Default configuration is RS-232 with pin 6 disconnected. Connect to pin 2 (0 V) to enable RS-485 serial communications.
7	Fail – status output/ Seri- al-TX/RS-485 B+	-	Logic high when a fail/fault condition exists and the serial enable control input is inac- tive.
8	0 V control reference	-	0 V reference for all control and status sig- nals listed within this table.
9	Analogue speed – con- trol input	-	8.4-10 V analogue input: 8.4 V = 83% speed; +10 V = 100% speed
10	Chassis/Screen	-	Screen
11	+10 V analogue refer- ence – Control output	Positive	+10 V analogue voltage reference output: 5 mA; unipolar output, diode protected.
12	Chassis/Screen	-	Screen
13	Not connected	-	Unused control pin
14	Remote – control input	-	Connect to pin 2 (0 V) to enable remote control via parallel or serial control modes.
15	Normal – status output	-	Logic low when the pump rotational speed is at normal speed or above.

6.4. Parallel control and monitoring



CAUTION: EMF RISK

Risk of damage to equipment. If you use the normal and fail lines to drive the coils of d.c. relays, include a back EMF suppression diode in parallel with each relay coil to protect the pump.

Connect the control equipment to the control input pins of the logic interface mating half. Refer to *Table: Logic interface pins* to identify the logic interface connector pins. The control inputs are:


- Start
- Standby speed
- Analogue speed
- A. To activate the control inputs, connect the relevant control input (pin 14) to the 0 V control reference.
- B. The NORMAL output can go down up to 100 mA when you control an external relay coil with an external coil voltage of + 24 V d.c. The external + 24 V voltage source must be referenced to the common control voltage of the pump control system, i.e. pin 2 of the 15-Way D-Type customer interface connector. Alternatively, if the NORMAL output is connected to + 10 V reference output of the control system, i.e. pin 11 of the 15-Way D-



Type connector, a 4.7 kOhm pull-up resistor is recommended to be kept in the current rating of + 10 V reference rail.

C. The FAIL output can go down up to 100 mA when you control an external relay coil with an external coil voltage of + 24 V d.c. The external + 24 V voltage source must be referenced to the common control voltage of the pump control system, i.e. pin 2 of the 15-Way D-Type customer interface connector. Alternatively, if the FAIL output is connected to + 10 V reference output of the control system, i.e. pin 11 of the 15-Way D-Type connector, a 4.7 kOhm pull-up resistor is recommended to be kept in the current rating of + 10 V reference rail.

Figure 7 Logic interface connections - parallel control



1. Start switch

- 3. Pump logic interface
- 2. Standby switch (optional)

6.5. Analogue speed control

The analogue speed input is a process control source which enables the pump to run at variable operating speeds. This speed control source is an alternative to standby speed control.






Figure 8 Logic interface connections - analogue speed control

- (125 Hz to 150 Hz)
- 3. Active speed range 83.00% to 100.00%
- 4. Active Voltage range 8.4 V to 10.00 V

Note:

0.1 V = 1% of default run speed. Voltages below 5 V will result in a clamped speed of 50% of full speed.



6.6. Hardware configuration

Use the 15-way D-type connector (*Figure: Pump features*, (6)) and configure the signal to enable the analogue speed control source:

- Connect the analogue speed enable control input (pin 1) to the 0 V control reference (pin 2).
- Connect a suitably calibrated analogue voltage source (0 to +10 V), for example, (DAC) to the analogue speed control input (pin 9).
- As an alternative, connect the output of a potentiometer referenced to the pump reference voltage (pin 11) to the analogue speed control input (pin 9). Refer to *Figure: Logic interface connections analogue speed control*. The 0 V rail of the external voltage source must be connected to the 0 V control reference (pin 2) of the pump controller.

6.7. Operational modes

- A +10 V input results in a mechanical running speed which is equal to 100% of the default run speed, that is 150 Hz.
- The minimum running speed provided by the analogue speed control source, is clamped at the minimum standby speed setting (approximately 83% of the default run speed).
- The maximum running speed provided by the analogue speed control source is clamped by the maximum standby speed setting (100% of the default run speed).

6.8. Start the pump



CAUTION: EXHAUST PIPELINE BLOCKAGE

Risk of damage to equipment. Do not operate the pump if the exhaust pipeline is restricted or blocked. The pump will not operate correctly and can be damaged.

- 1. Make sure that the vacuum system isolation valve is closed (if installed).
- 2. Make sure that the mains supply to the pump is isolated.
- 3. To connect the pump to the power supply use the power cord which is part of the delivery scope.
- 4. Make sure that the exhaust extraction system is not restricted and that the valves in the exhaust extraction system are open.
- 5. Supply the power.
- 6. Start the pump with the applicable control source that follows:
 - manual control mode push the start button, refer to Figure: Control panel, item 2
 - · parallel control mode use the start, stop control input
 - serial control mode use the start command
- 7. Open the vacuum system isolation valve (if installed).

6.9. Shut down the pump



WARNING: PUMP SUPPLY ISOLATION

Risk of injury. Do not remove the inlet connections until the pump stops rotating and the power and nitrogen supply has been isolated. The pump can take up to 3 minutes to fully stop.





CAUTION: ELECTRICAL SUPPLY

Risk of damage to equipment. Do not disconnect the pump from the electrical supply until the pump has fully stopped.

You can shut down the pump with:

- the front panel controls
- the 15-way D-type connector in either parallel or serial control mode.

Note:

If the pump is to be shut down for storage, remove any process gases by running on a gas ballast for at least one hour.

To shut down the pump, obey the instructions that follow:

- 1. Close the vacuum system isolation valves (where installed) to prevent the flow of the liquid back into the vacuum system.
- 2. Stop the pump system with the applicable control source that follows:
 - manual control mode push and hold the stop button (Figure: Control panel, item 3)
 - · parallel control interface use the start, stop control input
 - serial control mode send a stop command
- 3. When the pump is set to off, the fan will operate until the pump cools down.
- 4. Vent the pump through the inlet.
- 5. Isolate the mains supply.

6.10. Restart the pump

If the pump is shut down automatically because of high pump power make sure the pump rotates freely. Restart the pump. If the problem occurs again, contact us..

If the pump is shut down automatically because of an alarm condition, correct the alarm condition before you start the pump. To restart the pump, refer to *Start the pump* on page 30.

6.11. Unplanned shutdown

Allow the pump to cool after a sudden power outage, for example, due to power cut, circuit breaker trip or if the power cord is removed accidentally. We recommend a minimum of one minute to cool the pump. More time is required if the pump is working with a high load or high ambient temperature.


If the pump is configured for auto-restart, the pump may trip due to overcurrent if the power cut is less than 7 seconds. A longer cooling period is required to reset the pump. Disconnect the power supply for one minute and then reconnect it. If an error occurs, refer to *Fault finding*.

6.12. Gas ballast operation

If water vapor is present within the system, or will be pumped by the product, the gas ballast must be used to prevent damage to the product. The gas ballast should continue to be applied with the pump running, for at least 1 hour after the water vapor source has been removed to enable the pump to dry.



Figure 10 Maximum vapour handling capacity



6.13. Gas ballast connection



WARNING: HAZARDOUS SUBSTANCES


Risk of explosion. The pump is intended for clean applications only. It is not suitable for use with hazardous production materials.


Gas ballast is operated by Let-Lok pipe fitting. It enables externally actuated valves or a controlled gas supply to be used through a metal pipe with 8 mm outer diameter.

To install the pipe:

Remove the plug (item 1,2 and 3). Assemble nut, front and back ferrule. Tighten the nut by hand. Insert the pipe and check if it rests firmly on the fitting shoulder. Use a wrench to tighten the nut. One and quarter turn is required to assure an effective seal. Refer to *Gas ballast operation*.







7. Maintenance and service

7.1. General maintenance

The system requires no user maintenance. Maintenance must be done by our service centres. The frequency of the maintenance depends on the process. The frequency of the maintenance can be adjusted based on pump use.

We recommend a monthly visual inspection of the dry pump system. Check that the electrical supply cord, hoses and pipelines connected to the pump are in good condition and tightened.

7.2. Maintenance plan

More frequent maintenance may be required if the pump is used to pump aggressive gases or vapours, solvents, organic substances and acids, or if the pump is operated continuously at high operating temperature.

Operation	Frequency (months)	Service indi- cator	Reference
Inspect and clean the inlet strainer	12	No	Inspect and clean the inlet strainer
Inspect and clean the external fan cover if required	12	No	Clean the exter- nal fan grill
Electrical safety check	60	No	Electrical safety check

Table 13 Maintenance plan

7.3. Inspect and clean the inlet strainer

Whenever the pump is disconnected from the vacuum system or annually:

- Remove the inlet strainer from the pump inlet. Refer to *Figure: Pump features*.
- Clean the pump inlet and remove the debris that may have accumulated.
- Examine the inlet strainer. If necessary, clean it with a cleaning solution suitable for the substances pumped.
- Install the inlet strainer before you reconnect the pump to the vacuum system. Refer to *Mechanical installation*.

7.4. Clean the ventilation slots

If the ventilation slots are not kept clean, the air flow over the pump can be restricted and the pump may overheat.

To clean the ventilation slots:

- 1. Switch off the pump and disconnect from the electrical supply.
- 2. Use a dry cloth and a soft brush to remove dirt and deposits on the ventilation slots.

7.5. Electrical safety check

Test the earth continuity and the insulation resistance of the pump system in accordance with local regulations for the periodic test of electrical equipment.



The earth continuity must be less than 0.1 Ω and the DC insulation resistance greater than 1.0 M Ω . If the pump fails any of these tests, contact an field service employee (FSE).

7.6. Return the equipment or components for service (SCIEX)

When returning the pump for service or for any other reason, you must complete the required forms provided to you for Declaration of Contamination of Vacuum Equipment and Components. The competed form tells us if any substances found in the equipment are hazardous, which is important for the safety of our employees and all other people involved in the service of your equipment.



8. Fault finding

A list of fault conditions and their possible causes is given to assist in basic troubleshooting. If you are not able to correct a fault, call your supplier or your nearest manufacturer service centre for advice.

Fault The pump has failed to start or has stopped

Cause The electrical supply fuse has blown.

Remedy Make sure that the external electrical supply is set to on and that the power indicator LED (*Figure: Control panel*, item 14) is illuminated. If the LED does not illuminate, contact us.

Cause The motor controller has not reset after an over temperature event.

Remedy Disconnect the main power supply and wait for at least one minute and attempt to restart the pump. If the pump still does not start, contact us.

If the alarm indicator LED (*Figure: Control panel*, item 5) flashes, it indicates that the pump has high current load from possibly a high inlet pressure/flow. Adjust the inlet pressure/flow until the pump performance recovers.

Fault The pump has failed to achieve the required performance

- Cause The pressure measurement technique or gauge head is not suitable or gives an incorrect indication of pressure.
- Remedy Make sure that the vacuum measurement equipment is calibrated and is updated. The correct gauge range must be selected for the application. Contact us for further assistance on gauge selection, if required.
- Cause The vacuum fittings are dirty or damaged.
- Remedy Make sure that the vacuum fittings are clean and scratch free.
- Cause There is a blockage or high pressure in the exhaust line.
- Remedy Make sure that the exhaust valves fitted are not closed when the pump is in operation.

Cause The motor controller is current limiting the supply.

Remedy If the alarm indicator LED (*Figure: Control panel*, item 5) flashes, the pump has high current load, possibly from a high inlet pressure/flow. Adjust the inlet pressure/flow until the pump performance recovers.

Fault The pump is noisy

Cause The pump is contaminated with solid particles.

Remedy Contact us or the supplier for further information.

Fault The pumping speed is poor or pump down time is too long

Cause The pipelines connections are too small in diameter.

Remedy Make sure that the pipework has sufficient conductance (user's responsibility) and the pump performance is not compromised.

Cause There is a leak in the system.

Remedy Do the leak test of the pump system in accordance with the requirements specified in *Table: Performance data*.



Cause The motor controller is current limiting the supply.

Remedy Adjust the inlet pressure/flow until the pump performance recovers.

Cause The pump chamber is too big.

Remedy Make sure that the chamber size is not more than the limits specified in *Table: Mechanical data*.

8.1. Alarm indicator codes

When the fail condition becomes active, the red alarm indicator shows a flashing sequence.

- If the error light is on continuously, this indicates a problem has been found with the embedded software. In this case, switch the power supply on and off.
- If the indication is not cleared, a software download may be required. In this case, contact us or the supplier.
- If the alarm indicator is flashing, identify the error flash code. Refer to *Table: Flashing error codes*.

Note:

There is a sufficient off period between each subsequent cycle repetition to mark the start of a new flash sequence. The duration of a long flash (L) is equal to 3 times the duration of a short flash(s) (0.5 s).

Error flash po- sition	Error flash sequence	Comment	Actions
0	S-S-S-S-S-S	Overload timeout	Check if the pump is not under constant high pressure or the inlet or outlet is not blocked.
1	L-s-s-s-s	Controller software er- ror	Switch the power to the pump off and on and see whether the error code appears again. If it does, contact us or the supplier.
2	s-L-s-s-s-s	Controller failed inter- nal configuration and calibration operation	Switch the power to the pump off and on and see whether the error code appears again. If it does, contact us or the supplier.
3	s-s-L-s-s-s	Acceleration timeout	Check if the pump is under constant high pressure or the inlet or outlet is blocked.
4	s-s-s-L-s-s	Over current trip acti- vated, or other hard- ware fault	Switch the power to the pump off and on and see whether the error code appears again. If it does, contact us or the supplier.
5	s-s-s-s-L-s	Self test fault	Switch the power to the pump off and on and see whether the error code appears again. If it does, contact us or the supplier.
6	s-s-s-s-s-L	Serial control mode interlock	Reactivate the serial enable and send a se- rial command to clear the error code.

Table 14 Flashing error codes



9. Storage



WARNING: INHALATION HAZARD

Risk of asphyxiation. Do not burn the fluoroelastomer seals and O-rings.



CAUTION: LUBRICATION TRANSFER

Risk of damage to equipment. Install blanking plates to seal the vacuum inlet and outlet port. The pumps must be stored in a horizontal configuration to prevent possible lubrication transfer.

To store the pump:

- 1. Shut down the pump. Refer to *Shut down the pump* on page 30.
- 2. Disconnect the pump from the electrical supply.
- 3. Place and secure protective covers over the inlet and outlet ports.
- 4. Store the pump in a clean and dry condition until required for use.

When necessary, prepare and install the pump. Refer to *Installation* on page 16.

10. Disposal



WARNING: CONTAMINATION HAZARD

Risk of toxic exposure and acid burns. Identify, contain and safely dispose of contaminated items.

Dispose of the pump and any components or accessories safely and in accordance with all local and national safety and environmental requirements.

Dry pump system materials suitable for recycling include cast iron, steel, PTFE, stainless steel, aluminium, zinc alloy, nickel, mild steel, ABS and polyamide.

Take particular care with:

- Fluoroelastomers which may have decomposed as the result of being subjected to high temperatures.
- Components which have been contaminated with dangerous process substances.