

# Amino Acid Analysis Quick Reference Card

## For Hydrolysate Samples

### Safety

For safety and biohazard guidelines, refer to the “Safety” section in the Amino Acid Analysis for Hydrolysate Samples Protocol. For all chemicals in **bold red** type, read the MSDS and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves.

### Overview

This document supports amino acid analysis for hydrolysate samples using iTRAQ™ Reagents and the AB SCIEX Amino Acid 20/20 Analyzer.

The labeling protocol labels a peptide hydrolysate, protein hydrolysate, or a hydrolysate from animal feed sample (dry  $\approx 10$  nmol amino acid) with iTRAQ Reagent 117. An iTRAQ™ Reagent 114-labeled amino acid standard is added as an internal standard.

See the Amino Acid Analysis for Hydrolysate Samples Protocol for:

- Pipetting recommendations.
- Recommended analysis conditions, quality control, and troubleshooting.
- Supplemental information such as “How to Obtain More Information.”

### Amino Acid Analysis of Hydrolysate Samples - iTRAQ™ Reagents Labeling Protocol

**IMPORTANT!** Throughout the procedure, cap each tube promptly to avoid evaporation.

#### Labeling the Hydrolysate Sample with iTRAQ Reagent 117

1. If necessary, dry the hydrolysate sample.  
**IMPORTANT!** For optimal labeling, the hydrolysate sample must be completely dry.
2. Add Labeling Buffer - Amino Acid†:
  - To each sample tube containing 1  $\mu\text{g}$  hydrolysate:
    - a. Add 5  $\mu\text{L}$  Labeling Buffer - Amino Acid.
    - b. Vortex to mix, then spin.
  - To each sample tube containing more than 1  $\mu\text{g}$  hydrolysate:
    - a. For every 1  $\mu\text{g}$  of hydrolysate, add 5  $\mu\text{L}$  of Labeling Buffer - Amino Acid. For example, if your sample contains 6  $\mu\text{g}$  of hydrolysate, add 30  $\mu\text{L}$  of Labeling Buffer - Amino Acid
    - b. Vortex to mix, then spin.
3. Spin each required vial of **iTRAQ Reagent 117** (at room temperature) to bring the solution to the bottom of the vial.
4. Add 70  $\mu\text{L}$  of **isopropanol**. Mark the vial as “diluted.”

5. Vortex each vial to mix, then spin.
6. Label each hydrolysate sample prepared in [step 2](#) with iTRAQ Reagent 117:
  - To each sample tube containing 1  $\mu\text{g}$  hydrolysate:
    - a. Add 5  $\mu\text{L}$  of diluted iTRAQ Reagent 117.
    - b. Vortex to mix, then spin.
  - To each sample tube containing more than 1  $\mu\text{g}$  hydrolysate:
    - a. Transfer a 5- $\mu\text{L}$  aliquot of the hydrolysate sample/ Labeling Buffer - Amino Acid solution to a fresh tube.
    - b. To the aliquot, add 5  $\mu\text{L}$  of diluted iTRAQ Reagent 117.
    - c. Vortex to mix, then spin.
7. Incubate the sample tubes at room temperature for at least 30 min.
8. Add 1  $\mu\text{L}$  of hydroxylamine to each sample tube.
9. Vortex each sample tube to mix, then spin.
10. Incubate the sample tubes at room temperature for at least 5 min.
11. Dry the samples completely in a centrifugal vacuum concentrator (generally not more than 1 hour).

#### Adding iTRAQ Reagent 114-Labeled Internal Standard

1. Prepare a 6-pmol/ $\mu\text{L}$  iTRAQ™ Reagent 114-labeled amino acid internal standard by reconstituting one vial of **Hydrolysates Standards - 114 Labeled** with Sample Diluent - Amino Acid. The amount of Sample Diluent - Amino Acid to use is indicated on the Certificate of Analysis (approximately 1.67 mL).
2. Vortex to mix, then spin.  
The iTRAQ™ Reagent 114-labeled amino acid internal standard can be stored at  $-15$  to  $-25$  °C.
3. Add 25  $\mu\text{L}$  of the iTRAQ™ Reagent 114-labeled amino acid internal standard to each dried iTRAQ Reagent 117-labeled sample.
4. Vortex each tube to mix, then spin.

**IMPORTANT!** This procedure yields enough material for approximately three 5- $\mu\text{L}$  injections for each sample. Discard any remaining material.

### iTRAQ™ Reagent-Labeled Amino Acids in a 5- $\mu\text{L}$ Injection

A 5- $\mu\text{L}$  injection prepared according to the protocol contains:

- iTRAQ Reagent 117-labeled amino acids in the sample
- 30 pmol of iTRAQ Reagent 117-labeled norvaline
- 30 pmol of each iTRAQ Reagent 114-labeled amino acid in the standard, including norvaline

† Labeling Buffer - Amino Acid contains 30 pmol/ $\mu\text{L}$  norvaline.

## Equilibrating the Column Before Reuse

1. Wash the column with 2.5 mL 30% Mobile Phase A/ 70% Mobile Phase B at 1.0 mL/min for 2.5 min.
2. Equilibrate the column with 25 mL of 98% Mobile Phase A/ 2% Mobile Phase B (initial starting conditions) at 1.0 mL/min for 25 min.
3. To verify a stable baseline, continue the 98% Mobile Phase A/2% Mobile Phase B flow and inject 5  $\mu$ L Sample Diluent - Amino Acid (3.5% formic acid) at least twice.

## Mobile Phases, Gradient, and Chromatogram


**Table 1 Recommended mobile phase preparation (1L)**

Component	Mobile Phase A	Mobile Phase B
Mobile Phase Modifier A	1.00 mL	1.00 mL
Mobile Phase Modifier B	50.0 $\mu$ L	50.0 $\mu$ L
Milli-Q® water HPLC grade	998.95 mL	—
Acetonitrile, HPLC-grade	—	998.95 mL

**Table 2 Recommended LC gradient**

Total Time (min)‡	%Mobile Phase A	%Mobile Phase B
0.0	98.0	2.0
10.0	72.0	28.0
10.1	0.0	100.0
16.0	0.0	100.0
16.1	98.0	2.0
25.0	98.0	2.0

‡ Recommended flow rate: 1.0 mL/min, split in the source to 200 to 250  $\mu$ L/min



**Figure 1 Representative chromatogram**

## Amino Acid Specifications

For each period in the MRM experiment:

- Q1 for each iTRAQ Reagent 114-labeled amino acid in the internal standard and the corresponding iTRAQ Reagent 117-labeled amino acid in the sample is the labeled mass shown in the table.
- Q3 for each labeled amino acid in the internal standard = 114.
- Q3 for each labeled amino acid in the sample = 117.

Name	Formula	MH+ (amu)	
		Unlabeled	Labeled (Q1)
Period 1, Experiment 1 Dwell = 50			
L-serine	C <sub>3</sub> H <sub>7</sub> NO <sub>3</sub>	106.1	250.2
glycine	C <sub>2</sub> H <sub>5</sub> NO <sub>2</sub>	76.0	220.1
L-histidine	C <sub>6</sub> H <sub>9</sub> N <sub>3</sub> O <sub>2</sub>	156.1	300.2
L-aspartic acid	C <sub>4</sub> H <sub>7</sub> NO <sub>4</sub>	134.0	278.2
L-methionine sulfoxide	C <sub>5</sub> H <sub>11</sub> NO <sub>3</sub> S	166.1	310.2
L-threonine	C <sub>4</sub> H <sub>9</sub> NO <sub>3</sub>	120.1	264.2
L-alanine	C <sub>3</sub> H <sub>7</sub> NO <sub>2</sub>	90.1	234.2
L-glutamic acid	C <sub>5</sub> H <sub>9</sub> NO <sub>4</sub>	148.1	292.2
L-arginine	C <sub>6</sub> H <sub>14</sub> N <sub>4</sub> O <sub>2</sub>	175.1	319.2
Period 2, Experiment 1 Dwell = 100			
L-proline	C <sub>5</sub> H <sub>9</sub> NO <sub>2</sub>	116.1	260.2
L-cysteine	C <sub>3</sub> H <sub>7</sub> NO <sub>2</sub> S	122.0	266.1
L-lysine (2 labels)	C <sub>6</sub> H <sub>14</sub> N <sub>2</sub> O <sub>2</sub>	147.1	435.3
Period 3, Experiment 1 Dwell = 50			
L-valine	C <sub>5</sub> H <sub>11</sub> NO <sub>2</sub>	118.1	262.2
L-norvaline	C <sub>5</sub> H <sub>11</sub> NO <sub>2</sub>	118.1	262.2
L-methionine	C <sub>5</sub> H <sub>11</sub> NO <sub>2</sub> S	150.1	294.2
L-tyrosine	C <sub>9</sub> H <sub>11</sub> NO <sub>3</sub>	182.1	326.2
L-isoleucine	C <sub>6</sub> H <sub>13</sub> NO <sub>2</sub>	132.1	276.2
L-leucine	C <sub>6</sub> H <sub>13</sub> NO <sub>2</sub>	132.1	276.2
L-norleucine	C <sub>6</sub> H <sub>13</sub> NO <sub>2</sub>	132.1	276.2
L-phenylalanine	C <sub>9</sub> H <sub>11</sub> NO <sub>2</sub>	166.1	310.2

© 2010 AB SCIEX. The trademarks mentioned herein are the property of AB Sciex Pte. Ltd. or their respective owners. AB SCIEX™ is being used under license.

For Research Use Only. Not for use in diagnostic procedures.

08/2010