

Top/Middle-Down Protein Sequencing Demo Workflows

Biologics Explorer Software Guidelines

Powered by Genedata Expressionist[®]



TPUB-CUST-16032-A

© 2024 DH Tech. Dev. Pte. Ltd.



Top/Middle-Down Protein Sequencing Demo Workflows

CONTENTS OF THIS GUIDE

A: Overview of the Protein Sequencing Demo Workflows

1. Overview of the Applications for Protein Sequencing Demo Workflows
2. Overview of the Protein Sequencing Demo Workflows

B: Information about the Protein Sequencing Demo Workflows

1. Protein Sequencing
2. Protein Sequencing Review

Part A

Overview of the Protein Sequencing Demo Workflows

A

Overview of the Applications for Protein Sequencing Demo Workflows

- The Protein Sequencing demo workflows contain examples of how to analyze top-down or middle-down protein sequencing data of:
 - Antibody subunits. For example, Fd', Fc/2, and LC after IdeS digestion and reduction of disulfide bonds.
 - Antibody chains. For example, HC and LC after reduction of disulfide bonds.
 - Proteins with masses to a maximum of approximately 50 kDa.
- Data quality for protein sequencing can be increased if multiple consistent technical replicates are analyzed at the same time.
 - If the overall signal intensity is good, with a satisfactory signal-to-noise ratio, then a single sample can be sufficient for analysis.
- The **Top_Middle-Down_ProteinSequencing_Demo** workflow analyzes each antibody subunit or chain separately.
- The **Top_Middle-Down_ReviewSnapshots_Demo** reviews all of the results together.

Overview of the Protein Sequencing Demo Workflows

Top_Middle-Down_ProteinSequencing_Demo:

- A workflow for top-down characterization of small proteins or middle-down characterization of reduced and IdeS digested biotherapeutic molecules.

Top_Middle-Down_ReviewSnapshots_Demo :

- A workflow to review results of all subunits and show total sequence coverage from different fragmentation types.

Part B

Information About Protein Sequencing Demo Workflows


B

Protein Sequencing
Demo Workflow Information

B1

Overview and Application: Top_Middle-Down_ProteinSequencing_Demo


- This workflow uses data from 10 replicate injections of a biotherapeutic molecule after reduction and IdeS digestion.
- Each subunit (LC, Fd', Fc/2) is analyzed separately.
- MS/MS fragmentation data from the 3 selected charge states is summed together, and then averaged for each subunit.
- The search parameters in the *Fragment Mapping* activity node are optimized to identify MS/MS fragments, with post-translational modifications, and glycosylation at fixed positions.
- For more information about how to use Biologics Explorer software, refer to the *Biologics Explorer Quick Guide* and *Top/Middle-Down Protein Sequencing Template Workflow Guidelines*.



Protein Sequencing Review
Demo Workflow Information


B2

Overview and Application: Top_Middle-Down_ReviewSnapshots_Demo




- This workflow uses Snapshot files of results from middle-down analysis of a biotherapeutic molecule after reduction and IdeS digestion. Each subunit was fragmented with EAD and CID across 3 charge states.
- The sequence coverage for each subunit (LC, Fd', Fc/2) is the combination of EAD and CID fragment ions.
- For more information about how to use Biologics Explorer software, refer to the *Biologics Explorer Quick Guide* and *Top/Middle-Down Protein Sequencing Template Workflow Guidelines*.

Review Results



- **Combine by Proteoform** shows the combined sequence coverage of data acquired with EAD and CID fragmentation for each subunit (LC, Fd', Fc/2).

- Use the **Wrap** icon to compare Coverage for different sequences
- Use the **Toggle spectrum direction** icon to create a mirror plot in the **Fragment Spectra Viewer**.
- To change **Labels**, right-click the **Fragment Spectra Viewer**, and then select **Settings**.
 - To increase the number of labeled peaks that are shown, select **Label Features: Offset**.
- Annotations that are selected in **Coverage**, **Fragment Spectra Viewer** and **Fragment Table** are synchronized.



Trademarks/ Licensing

The SCIEX clinical diagnostic portfolio is For In Vitro Diagnostic Use. Rx Only. Product(s) not available in all countries. For information on availability, please contact your local sales representative or refer to www.sciex.com/diagnostics. All other products are For Research Use Only. Not for use in Diagnostic Procedures.

Trademarks and/or registered trademarks mentioned herein, including associated logos, are the property of AB Sciex Pte. Ltd. or their respective owners in the United States and/or certain other countries (see www.sciex.com/trademarks).

© 2024 DH Tech. Dev. Pte. Ltd. TPUB-CUST-16032-A





The power of precision

The SCIEX clinical diagnostic portfolio is For In Vitro Diagnostic Use. Rx Only. Product(s) not available in all countries. For information on availability, please contact your local sales representative or refer to www.sciex.com/diagnostics. All other products are For Research Use Only. Not for use in Diagnostic Procedures.

Trademarks and/or registered trademarks mentioned herein, including associated logos, are the property of AB Sciex Pte. Ltd. or their respective owners in the United States and/or certain other countries (see www.sciex.com/trademarks).

© 2024 DH Tech. Dev. Pte. Ltd. TPUB-CUST-16032-A

