



---

# Capillary Electrophoresis System

## Installation Qualification 1

---

This document is provided to customers who have purchased SCIEX equipment to use in the operation of such SCIEX equipment. This document is copyright protected and any reproduction of this document or any part of this document is strictly prohibited, except as SCIEX may authorize in writing.

Software that may be described in this document is furnished under a license agreement. It is against the law to copy, modify, or distribute the software on any medium, except as specifically allowed in the license agreement. Furthermore, the license agreement may prohibit the software from being disassembled, reverse engineered, or decompiled for any purpose. Warranties are as stated therein.

Portions of this document may make reference to other manufacturers and/or their products, which may contain parts whose names are registered as trademarks and/or function as trademarks of their respective owners. Any such use is intended only to designate those manufacturers' products as supplied by SCIEX for incorporation into its equipment and does not imply any right and/or license to use or permit others to use such manufacturers' and/or their product names as trademarks.

SCIEX warranties are limited to those express warranties provided at the time of sale or license of its products and are SCIEX's sole and exclusive representations, warranties, and obligations. SCIEX makes no other warranty of any kind whatsoever, expressed or implied, including without limitation, warranties of merchantability or fitness for a particular purpose, whether arising from a statute or otherwise in law or from a course of dealing or usage of trade, all of which are expressly disclaimed, and assumes no responsibility or contingent liability, including indirect or consequential damages, for any use by the purchaser or for any adverse circumstances arising therefrom.

**For research use only.** Not for use in diagnostic procedures.

AB Sciex is doing business as SCIEX.

The trademarks mentioned herein are the property of AB Sciex Pte. Ltd. or their respective owners.

AB SCIEX™ is being used under license.

© 2017 AB Sciex



AB Sciex LLC  
500 Old Connecticut Path  
Framingham, Massachusetts 01701  
USA

# Contents

---

<b>Chapter 1 Installation Preparation</b> .....	<b>4</b>
Customer Information.....	4
System Information.....	4
Pre-Installation Requirements.....	5
<b>Chapter 2 Installation Qualification</b> .....	<b>7</b>
Installation Qualification Checklist.....	7
<b>Chapter 3 Installation Qualification Data Log</b> .....	<b>9</b>
Corrective Actions.....	9
<b>Chapter 4 Signoff</b> .....	<b>11</b>
<b>Appendix A Topics for Familiarization</b> .....	<b>12</b>
<b>Revision History</b> .....	<b>13</b>

# Installation Preparation

# 1

---

**Note:** For regulatory and safety information for the capillary electrophoresis system, refer to the *Safety Notices*, *System Overview*, or *Operator Guide*.

---

**Note:** Before starting the installation, verify that the customer has completed the checklist in the *Site Planning Guide* or, for some systems, the *Preinstallation Manual*.

---

## Customer Information

Organization			
Address			
City			
Country			
Telephone		Zip/Postal code	
Site planner contact name		Fax	
E-mail address			

## System Information

Table 1-1 System Information

Model	
Location	
Serial number	

**Table 1-2 Detector Information**

(If installed) UV detector serial number	
(If installed) PDA detector serial number	
(If installed) LIF detector serial number	
(If installed) LIF laser serial number	

**Table 1-3 Controller Information**

Model	
Serial number	
Operating system	

**Table 1-4 Software Information**

32 Karat software version	
32 Karat software license serial number	
Firmware version	

## Pre-Installation Requirements

Before starting the installation, refer to the *Site Planning Guide* or, for some systems, the *Preinstallation Manual*, and verify that the customer has prepared the site for installation.

In [Table 1-5](#), list any requirements that have not been met.

**Table 1-5 Nonconforming Requirements**

Requirement	Discrepancy	Resolution/Comments

## Installation Preparation

---

**Table 1-5 Nonconforming Requirements (continued)**

Requirement	Discrepancy	Resolution/Comments

# Installation Qualification

# 2

Installation Qualification (IQ) confirms that the system and its modules and accessories have been supplied as ordered and that the system is properly installed.

The following table lists the steps to complete a general installation qualification. Follow the steps and then record the results below and in the [Installation Qualification Data Log on page 9](#). Record any corrective actions in [Corrective Actions on page 9](#).

## Installation Qualification Checklist

Task	Complete
<b>Order Verification</b>	
<ol style="list-style-type: none"><li>1. Obtain the purchase order or shipping document and then record the numbers in the <a href="#">Installation Qualification Data Log on page 9</a></li><li>2. Make sure that all of the ordered items are present and undamaged. If items are missing or damaged, then follow standard SCIEX procedures.</li></ol>	
<b>Instrument Installation</b>	
<ol style="list-style-type: none"><li>1. Unpack and install the system using the procedure in the <i>Service Guide</i>.</li><li>2. Verify the fuse setting for the voltage provided and then install the power cords.</li></ol>	
<b>Instrument Initialization</b>	
<p>Turn on the system and then verify that all systems, modules, and accessories respond properly, including passing all power on self-tests (POST).</p> <p>If the system fails any of the self-tests, then turn off the power and restart it. If the system fails the self-tests again, then notify the customer and implement corrective action. Describe the corrective action in <a href="#">Corrective Actions on page 9</a>.</p> <p>If the system still fails the self-tests, then terminate the installation qualification and notify the Territory Specialist.</p> <hr/> <p><b>CAUTION: Potential System Damage. Make sure that there are no sample trays present in the system. If sample trays are present on the rails during initialization, a collision might occur inside the system.</b></p> <hr/>	

## Installation Qualification

---

Task	Complete
<b>Safety Information</b>	
Inform the user of the location of the <i>Safety Notices</i> document and then note its location in the <a href="#">Installation Qualification Data Log on page 9</a> .	
<b>Note:</b> Some systems do not have a <i>Safety Notices</i> document. For those systems the safety information is in the <i>User Manual</i> , <i>System Overview</i> , or <i>Operator Guide</i> .	
<b>System Documentation</b>	
1. Inform the user of the location of all user documentation and then note its location in <a href="#">Installation Qualification Data Log on page 9</a> .	
2. Demonstrate how to access the software Help system.	
<b>Scheduled Maintenance</b>	
1. Inform the user that annual maintenance is recommended.	
2. Explain how to obtain emergency service and how to schedule planned maintenance. Note this information in the <a href="#">Installation Qualification Data Log on page 9</a>	
<b>Operator Instruction</b>	
Familiarize the user with the items in <a href="#">Topics for Familiarization on page 12</a> .	



# Installation Qualification Data Log

# 3

---

Purchase order number	
Location of <i>Safety Guide</i> or other safety information	
Location of user documentation	
FSE contact information	

## Corrective Actions

If no corrective actions were required, check here and do not complete the following.	
---	--

Test name	
Corrective action	
Test name	
Corrective action	
Test name	

## Installation Qualification Data Log

---

Corrective action	
Test name	
Corrective action	
Test name	
Corrective action	

# Signoff

# 4

---

Organization			
Customer contact name		Date (yyyy-mm-dd)	
Customer contact signature*			
FSE name		Date (yyyy-mm-dd)	
FSE signature			

\* Signature required on hard copy only.

# Topics for Familiarization

---

# A

Familiarize or review the following with the customer:

- Chemical and biological safety
- Electrical safety
- Windows log in
- Software functions:
  - License
  - System administration
  - Creating, editing, saving, and printing a method
  - Opening data files
  - Creating, editing, and saving a sequence
  - Running a single sample or a sequence
  - Stopping a run
  - Displaying data
  - Creating and viewing custom reports
- Installing a cartridge
- Loading samples
- Changing detectors
- Maintenance procedures

# Revision History

---

Revision	Reason for Change	Date
144867AA	First release of document.	September 1998
144867AB		January 1999
144867AC	Changed document title.	October 1999
144867AD		November 2000
144867AE		November 2011
RUO-IDV-08-5454-A   144867AF	Re-branded. Applied new template and converted to a PDF form.	May 2017