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Intact Analysis of Antibody Drug Conjugates

Trastuzumab Emtansine Analysis using Benchtop X500B QTOF Mass Spectrometer

Wen Jin¹, Suya Liu¹, Doug Simmons¹, Ian Moore¹, Sean McCarthy² and Sibylle Heidelberger¹

71 Four Valley Dr. Concord, ON L4K 4V8, Canada

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Column	Agilent Poroshell 300SB-C8 1.0 x 75mm 5 μm
Mobile phase A	0.1 % Formic acid in water
Mobile Phase B	0.1% Formic acid in acetonitrile
Flow rate	0.2 mL/min
Column Temperature	75 °C

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Source parameters:	
CUR	30
GS1	50
GS2	50
Ion Spray Voltage	5000 V
Source Temperature	400°C
TOFMS mass range	900 – 4000 m/z
DP	250 V
Accumulation time	0.5 s
Time bins to sum	80
Intact Protein Mode (IPM)	On
Large Proteins (>70kDa)	On
Decrease Detector Voltage	On
DP	250 V

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- [1] Michael T. Kim, Yan Chen, Joseph Marhoul and Fred Jacobson. Statistical Modeling of the Drug Load Distribution on Trastuzumab Emtansine (Kadcyla), a Lysine-Linked Antibody Drug Conjugate. *Bioconjugate Chem.* 2014, 25, 1223-1232
- [2] Yan Chen, Michael T. Kim, Laura Zheng, Galahad Deperalta, and Fred Jacobson. Structural Characterization of Cross-Linked Species in Trastuzumab Emtansine (Kadcyla). *Bioconjugate Chem.* 2016, 27, 2037-2047
- [3] Liuxi Chen, Lan Wang, Henry Shion, Chuanfei Yu, Ying Qing Yu, Lei Zhu, Meng Li, Weibin Chen, and Kai Gao. In-depth structural characterization of Kadcyla (ado-trastuzumab emtansine) and its biosimilar candidate. *MABS*, 2016, VOL. 8, NO. 7, 1210-1223

For Research Use Only. Not for use in Diagnostics Procedures.

AB Sciex is operating as SCIEX.

© 2019. AB Sciex. The trademarks mentioned herein are the property of AB Sciex Pte. Ltd. or their respective owners. AB SCIEX™ is being used under license.

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