ekspert ultraLC 100-XL

Release: 1.0 SP#1

Sep 2012

AB SCIEX Headquarters 500 Old Connecticut Path Framingham, MA 01701

U.S.A

Technical Support: 1 800-952-4716 or 877-740-2129

E-mail: Support@absciex.com

Fax: 1 650-627-2739

Internet: http://www.absciex.com

Copyright notices

ekspert ultraLC 100-XL

Copyright(c) 2012 AB SCIEX. All rights reserved.

What is in this document

Required firmware.

Redistribution and copies of the software.

Minimum PC configuration.

Scanned for viruses.

Uninstall notes.

Minimum required firmware:

ekspert 100-XL Autosampler: Flash 0950.110-1.02 (version 1.02)

Template 1 0013.115-01 Template 2 0013.116-01

ekspert 100 Pump: Flash 0760.100-1.12 (version 1.12) ekspert 100 Column Oven: Flash 0886.100-2.04 (version 2.04)

Redistribution and copies of the software:

You can make as many copies of the software as you wish. You are also allowed to redistribute copies of the software.

Minimum pc configuration:

ekspert ultraLC 100-XL runs on

Microsoft Window XP-Pro (SP2) operating systems and Microsoft Windows 7 operating systems with .Net framework 4.0 (client profile) installed. Furthermore, Analyst 1.6 must be installed. The following PC configuration is the minimum for ekspert ultraLC 100-XL:

- Pentium III 1.0 GHz (1000 MHz).
- 512 Mb internal memory.
- Hard disk with 850 MB of free storage space (less when .Net framework is installed already)
- SVGA video driver 1024 x 768; 16 bit "high color";
- Serial communication port for the system controlled.
- Mouse.
- CD-ROM drive.

It is strongly recommended to install the software on a local hard disk.

1024 Mb internal memory is strongly recommended for a PC that also runs other tasks like a virus scanner, remote access software and/or printer tasks.

Scanned for viruses:

The software is scanned for viruses by means of

- ESET NOD32 Antivirus version 4.2.40.0
- Virus signature database: 7450 (20120905)
- Update module: 1040 (20120313)
- Antivirus and antispyware scanner module: 1365 (20120724)
- Advanced heuristics module: 1134 (20120903)
- Archive support module: 1150 (20120828)
- Cleaner module: 1057 (20120626)
- Anti-Stealth support module: 1032 (20120806)
- ESET SysInspector module: 1214 (20100127)
- Self-defense support module: 1018 (20100812)
- Real-time file system protection module: 1006 (20110921)

Uninstall notes:

When uninstalling this software, only those files are removed that were installed.